

# CDF 200 - Prot. - Termofilme

|  |   |                   |                          |                        |                          |               |                   |               |                |
|--|---|-------------------|--------------------------|------------------------|--------------------------|---------------|-------------------|---------------|----------------|
| <b>File Name</b>                           | C:\USERS\MARCELO\DOCUMENTS\FTIR\CDF200P1.SP | <b>Date Stamp</b> | Mon Sep 16 14:24:06 2013 | <b>Date</b>            | Mon Sep 16 14:24:06 2013 |               |                   |               |                |
| <b>Technique</b>                           | Infrared                                    | <b>Instrument</b> | 16PC                     | <b>Spectral Region</b> | IR                       | <b>X Axis</b> | Wavenumber (cm-1) | <b>Y Axis</b> | %Transmittance |
| <b>Spectrum Range</b> 400.0000 - 4000.0000 |   |                   | <b>Points Count</b>      | 3601                   | <b>Data Spacing</b>      | 1.0000        |                   |               |                |

